



Personal incubators for separation of sensitive samples to reduce cross contamination risks

The new MiniGalaxy Personal CO<sub>2</sub> Incubator continues the fine manufacturing tradition established in the RS Biotech Galaxy series. Two models are available, the more sophisticated Mini A which supports a wide range of options, and the entry level Mini E.

small footprint with 50ltr capacity

Occupying only 48 x 47cm of bench space, the small footprint allows these incubators to be truly 'personal', thereby eliminating the risk of cross-contamination from cell lines of other users.

### ordering information

Mini CO <sub>2</sub> incubator 0.2-20% CO <sub>2</sub>	order code	order code
voltage 220/240V 50/60Hz	050-000	050-500
100/120V 50/60Hz	050-001	050-501
<b>No charge option (factory fitted)</b>		
Left Opening Door (hinged at right)	050-027	050-527
<b>Model options (factory fitted)</b>		
O <sub>2</sub> Control 0.1 - 19%	050-031	050-531
O <sub>2</sub> Control 1 - 19%	050-031	050-531
O <sub>2</sub> Control 1 - 95%	050-032	050-532
High Temperature Decontamination	050-034	050-534
Humidity Measurement & Display (includes all Mini A models)	050-037	050-537
Humidity Warning System	050-038	050-538
Access Port, diameter 25mm	050-019	050-519
RS 232 Port*	050-020R	050-520R
Centronics Port*	050-020C	050-520C
Soft Inner Door*	050-022	050-522
BMS Relay Contact Alarm*	050-023	050-523
Internal IP68 Sealed 240V Socket	050-024	050-524

\* Can be retrofitted by Service Personnel

### accessories

Galaxy IR CO <sub>2</sub> Analyser	order code
Galaxy Combined CO <sub>2</sub> and O <sub>2</sub> analyser	010-100
Perforated Shelf	010-1000
Non-perforated Shelf	050-104
Floor Stand, 400mm height for double stacking (part no. 050-118 also required)	050-105
Galaxylog V1 0.1°C based data collection software and RS232 port*	050-110
Stacking Frame for double stacking	050-117
Multi-position shelf racks - maximum 8 positions	050-118
Automatic CO <sub>2</sub> Change over Unit	050-121
CO <sub>2</sub> Two Stage Regulator	170-101
CO <sub>2</sub> Filter (supplied in pairs)	170-102
In-line Pressure Regulator	170-103
CO <sub>2</sub> Gas Analyser complete with 10 Testing Tubes	170-110
10 spare Gas Analyser Testing Tubes	170-113
Air Zero Filters (pack of 5)	170-114
	170-115

RS Biotech recommends that all in-line pressure regulators (part no. 170-110) should be set to gas supply from whatever O<sub>2</sub> source is employed

\* warnings with options a, Level 2 is



MiniGalaxy Incubators. Personal space at your fingertips

one space  
one user  
no sharing.....  
peace of mind

C&M SCIENTIFIC  
Unit 4, Lindsay Square, Deans Industrial Estate, Livingston, West Lothian, EH54 8RL  
Tel: 08452 608 601 Fax: 08452 608 602

