

LABORATORY FREEZER



Ask for a
Quotation

CM583WF

- Suitable for reagent, sample storage.
- Temperature range: -15/-25°C.
- Volume, gross: 583 ltr.
- 5 Year parts & labour warranty.

We're Scotland's leading scientific calibration, maintenance and validation company providing engineering support and consultation to a variety of organisations. Based in Livingston, we see your business as our business, helping our clients maintain their research platforms, adhere to accreditation standards and provide maintenance & support on temperature controlled equipment, centrifuges and more.

Accredited to ISO 9001, ISO 17025 & ISO 14001 as well as adhering to FDA and MHRA standards we consider ourselves your one stop shop for instrument servicing, repair and replacement of laboratory equipment.

Our team of engineers are trained to the highest level to deliver our wide range of technical services. We provide a call out service and welcome your enquiries by telephone or email. Get in touch to learn more about how we can service and support your laboratory operations. As well as how we can consult with you to address any issues or regulatory requirements you may have.

UK Head Office
1 Brewster Sq
Brucefield Industrial Estate
Livingston
West Lothian
EH54 9BJ

Tel: 01506 463 734
Email: info@cmscientific.co.uk
Web: www.cmscientific.co.uk

Overview

Suitable for reagent, sample storage. The 583 models are versatile freezers designed for storing products in laboratory workplaces.

Features

- Equipped with a right hand hinged reversible door with lock and 4 white shelves and castors.
- Units will meet all laboratory standards including ISO 15189 and cGMP when thermally mapped.
- Thermal mapping/validation service available on request.
- UKAS 17025 calibration available.
- Designed and built in Denmark for C&M Scientific Ltd.

Options

- Stacking kit for 125 Models.
- Castors on stacking kits.
- 30mm Porthole for sensor access.
- Drawer systems available on 125, 346, 583 models.

Specifications

| Variable | Value |
|---------------------------------|---|
| Model number | CM583WF |
| Category | Storage Cabinets |
| Type | Freezer |
| Line | Upright |
| Product Model | CM583WF |
| Equipped with | Right hand hinged reversible door with lock, 4 shelves. |
| Exterior | White |
| Interior | ABS lining in grey |
| Temperature range | -15/-25°C |
| Volume, gross | 583 ltr |
| Net usable volume | 436 ltr |
| Ex Dim. W x D x H, mm | 695 x 868 x 2000 |
| Weight, packed | 113 kg |
| Refrigerant / GWP | R 290 (propane)/3 |
| Ref.capacity at -10° | 512 Watt |
| Connection | 230 V, 50 Hz |
| Connection load | 485 Watt |
| Energy Efficiency | Class D |
| Energy consumption / year (AEC) | 2566 kWh |
| Sound level | 48.6 dB(A) |