

Ask for a
Quotation



CMG1350

CHROMATOGRAPHY CABINETS

- 600 Litre Cabinet.
- Extra Shelves & Bespoke Shelves.
- Bespoke Racking.
- Internal Sockets.
- Bespoke Internal Furniture.

We're Scotland's leading scientific calibration, maintenance and validation company providing engineering support and consultation to a variety of organisations. Based in Livingston, we see your business as our business, helping our clients maintain their research platforms, adhere to accreditation standards and provide maintenance & support on temperature controlled equipment, centrifuges and more.

Accredited to ISO 9001, ISO 17025 & ISO 14001 as well as adhering to FDA and MHRA standards we consider ourselves your one stop shop for instrument servicing, repair and replacement of laboratory equipment.

Our team of engineers are trained to the highest level to deliver our wide range of technical services. We provide a call out service and welcome your enquiries by telephone or email. Get in touch to learn more about how we can service and support your laboratory operations. As well as how we can consult with you to address any issues or regulatory requirements you may have.

UK Head Office

1 Brewster Sq
Brucefield Industrial Estate
Livingston
West Lothian
EH54 9BJ

Tel: 01506 463 734

Email: info@cmscientific.co.uk

Web: www.cmscientific.co.uk

Our chromatography cabinets are designed to accept chromatography instrumentation, and general purpose laboratory applications including sample storage in a temperature controlled environment.

Options

- 600 Litre Cabinet.
- Extra Shelves.
- Bespoke Shelves.
- Bespoke Racking.
- Internal Sockets.
- Bespoke Internal Furniture.

Specifications

Variable	Value
Model number	CMCG1350
Category	Storage Cabinet
Type	Chromatography cabinet
Line	Upright
Product Model	CMCG1350
Exterior / Interior	Stainless Steel & Aluminium
Temperature range	1°C - 10°C
Volume, gross	1350 ltr
Fan assisted cooling	2 x Internal Fans
Ex Dim. W x D x H, mm	1440 x 820 x 2080
Weight, packed	183 Kg
Refrigerant	134A
Fan assisted cooling	2 x Internal Fans
Shelves	6
Door	2 x Lockable spring loaded glass doors
Light	LED door lighting
Alarms	Hi / Lo Temp & Door open
Access Port	2 x 30mm
Internal sockets	2 x IP66 Double sockets